



HealthFriend IoT juicer

HOW TO USE HealthFriend IoT Juicer & app

MOTIV 1

*Children below age 7 and pregnant women are not allowed to use.
 *Not recommended for people with metal bone implants or pacemaker.

1

Click menu list, choose "Guest measurement".

2

Input basic information (age, gender, height, weight), then go to "Next".

3

After you read this content, go to "Next".

6

Firmly hold the upper and lower metal parts of the gauge at the same time.

5

Select measuring device then the measurement begins.

4

Press the (signal) button on juicer and check the LED light on, and select "Measure" button on app.

7

Complete measuring, ① drag up ② choose [Go to Print]. ③ select [Printing].

8

Check the print button and select it.

9

Take your result paper from printer.